


# Improved STR Artifact Recognition in Mixtures Using OSIRIS v. 2.8




George Riley, Robert Goor, Douglas Hoffman, Stephen Sherry  
National Center for Biotechnology Information, National Library of Medicine, NIH, Bethesda, MD

## Data analysis




Pull-up causes characteristic signals in raw data. Osiris analyzes these, along with cross-channel alignment, noise and other factors to determine the pattern of pull-up in the entire sample.

## OSIRIS Analyzes Partial Pull-up



OSIRIS calculates corrected RFU for alleles that contain some pull-up (15.2 is 1846 corr. RFU).

## OSIRIS Analyzes the Sample-wide Pull-up Pattern



Similar-sized peaks cause similar-sized pull-up in other channels. In this sample the pattern is shown with solid arrow indicating pull-up and dashed arrows and circles indicating lack of pull-up.


PU – pull-up, solid arrows; PPU – peak causing pull-up; Other artifacts not labeled. Identifier, analyzed data.

Pull-up is difficult to determine in mixtures with minor alleles based solely on single occurrence - whether a peak is truly pull-up or simply comigrates by coincidence. Pull-up shows a sample-wide pattern - peaks with equivalent amounts of DNA will cause equivalent pull-up.

OSIRIS v. 2.8 checks for a pull-up pattern and which are outlier peaks – alleles – in the pattern (least median of squares). OSIRIS removes alleles from pattern analysis, then uses regression to determine pull-up pattern. Off-scale data exhibit a different pattern and are considered separately.


This gives an accurate assessment of which peaks are pull-up (no allele call), and which are alleles (with pull-up correction for comigrating alleles), reducing editing and freeing up time for case interpretation.

## Pull-up Analysis in Mixtures




Automated pull-up and artifact determination facilitates efficient and accurate mixture analysis.

## More Osiris Improvements




Click to edit




Sample notices table

Click to edit alleles & artifacts in the Graph View. Change labels, add custom labels.

Scroll through and zoom to loci in the Sample notices table in the Graph View.



User-defined artifact labels and label priority



Artifact-specific labels

## OSIRIS Software

- Open source
- Validated as an expert system
- Cross-platform
  - PC (XP/Win7/8/10)
  - Mac OS10.8-10.10
- Displays
  - Stacked graph
  - Table
- Quality analysis
- Reanalysis prediction
- Flexible User configuration
- Flexible export
  - Table, LIMS, graphical
- Configured for most kits
- Files in both .hid and .fsa format
- Rapid analysis
  - < 30 s/96 samples

## Using OSIRIS

OSIRIS is a freely available download on the Osiris web page : <http://www.ncbi.nlm.nih.gov/projects/SNP/osiris/>

The Osiris User's Guide is on the homepage and in the program (F1) with a tutorial. The download includes demonstration data from various kits.

### Open Source Collaboration

OSIRIS source code is on the GitHub repository. We invite new collaborators to join the Osiris community to help improve Osiris.

For questions or requests, please contact: [forensics@ncbi.nlm.nih.gov](mailto:forensics@ncbi.nlm.nih.gov)

For announcements about Osiris releases and tips, subscribe to the Osiris announcement list from our homepage or: <http://www.ncbi.nlm.nih.gov/mailman/listinfo/forensics-announce>



Email a question.



Osiris Home page.

We would like to thank our collaborators at NIST and Jo Bright of ESR for their time discussing OSIRIS and data they provided. We thank Jalinnda Hull for her always excellent administrative assistance. Work at NCBI is supported by the Intramural Research Program of the NIH, National Library of Medicine.